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Manufacturer:
First Scientific Dental GmbH, 25335 Elmshorn, Germany

Esthetic low shrink light cure composite with nano-filler technology

EXCIFILL NANO is a **visible light cure nano-filled composite** for anterior and posterior restorations. Beside its low shrink and low abrasion EXCIFILL NANO shows improved handling features and excellent mechanical properties. EXCIFILL NANO is based on urethanedimethacrylate resin and inorganic filler particles <1.0 µm. The total filler content is 82 % (wt) and 74% (vol). The composite material meets the requirements of **DIN EN ISO 4049**.

Indications

Suitable for all cavities.

Side Effects

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of calciumhydroxide material.

Incompatibility with Other Materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the composite. Neither store the composite material in proximity of eugenol containing products, nor let the composite allow coming into contact with materials containing eugenol.

Preparation of Tooth and Colour Matching

Clean the tooth with a fluoride-free polishing paste (e.g. flour of pumice with water) prior to preparation and colour matching.

Ascertain the tooth shade while teeth are still moist and select the appropriate EXCIFILL NANO shade.

EXCIFILL NANO is shaded according to VITA®-Shades. It is also available in special shades and bleach shades.

Cavity Preparation

After isolation (best with a rubber dam) prepare the cavity with minimal tooth reduction.

For pulp protection areas close to the pulp should be covered with a thin layer of calciumhydroxide material.

Condition and apply a bonding agent according to the manufacturer instructions. The total etch technique with an etching gel and subsequent application of a light cure bonding agent (e.g. EXCIBOND) is recommended. Otherwise use a self etching bonding agent (e.g. EXCIBOND SE) before application of the EXCIFILL NANO.

Application and Curing of EXCIFILL NANO

Apply EXCIFILL NANO in the selected shade. EXCIFILL NANO is easy to model. Place it by using a suitable instrument. Matrix strips may be used.

In case of deep and extended cavities best results are obtained with a build-up of a base liner using a light cure resin modified cement or a light cure flowable composite, followed by the application of EXCIFILL NANO.

In case of small cavities EXCIFILL NANO can be applied directly.

For an optimum result apply EXCIFILL NANO in **layers of max. 2 mm**. Light cure each increment separately.

By using a polymerization unit (wavelength range 400-500 nm) with a **light intensity** of at least **1000mW/cm²**, cure each increment as follows:

Lighter shades (e.g. A1, A2, C2) **20 seconds**

Darker shades (e.g. A3.5, B3) **30 seconds**

Finishing

Depending on the shape of the area contour with finishing diamonds, flexible grinding discs or carbide burs.

Polish to high gloss with polishing discs or polishing brushes.

Storage

Do not store above 25°C (77°F). Protect from direct sunlight. Do not use after expiry date.

Close syringe immediately after use to avoid exposure to light.

Additional Notes

- Do not use any resin to adjust viscosity of composite restorative material.
- Contact of resin pastes with skin should be avoided, especially by anyone having known resin allergies.
- VITA® is a registered trade mark of the VITA-Zahnfabrik, Bad Säckingen, Germany.

Warranty

First Scientific Dental GmbH warrants this product will be free from defects in material and manufacture. First Scientific Dental GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and First Scientific Dental GmbH's sole obligation shall be repair or replacement of the First Scientific Dental GmbH product.

Limitation of Liability

Except where prohibited by law, First Scientific Dental GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Keep away from children!

For dental use only!

CE 0482

Manufacturer:

First Scientific Dental GmbH, 25335 Elmshorn, Germany